


MX4 Sort and Scan Instrument

Preliminary Product Information



MX4 Tube Management System

- Semi-automated mini-tube sorting/cherry-picking system
- Built-in full rack 2D barcode scanner
- Four source rack positions
- One destination rack position
- Work list based pick and place system
- Compatible with Matrix, Micronic, and ABgene 2D barcoded mini-tubes

XL Work List Manager Software

- Accepts work list of tube pick and place positions
- Guides operator through tray set up process
- Continually scans tubes for error checking
- Provides output file of destination rack content and tube position
- User defined error-handling
- LIMS system integration ready

www.biomicrolab.com

BioMicroLab, Inc.
2500-A Dean Lesher Drive
Concord, CA 94520
USA

Tel: 925-689-1200
Fax: 925-689-1263
Email: info@biomicrolab.com
Web: www.biomicrolab.com