

## Rack Thawing Station

### Defrosts Tube Racks in less than 20 minutes

Compatible with ABgene, MATRIX, & Micronic 96 tube rack labware



#### Unit Specifications:

- Dim: 6.25" W x 4.25" D x 11" H
- Weight: 5 lbs.
- Electrical: 110-220 VAC 50/60HZ

#### Rack Thawing Station

##### Fast performance

- Defrosts 96 tubes in < 20 minutes
- No warm up time

##### Small footprint

- Uses little bench top space

##### Best price

- Under \$599.00

