


Reagent Monitor SystemTM



Order PC based reagent level monitor & alarm

Applications

- Verify reagent start-up volumes
- Reduce low volume down-time issues
- SOP & QC task confirmation
- Monitor reagent use & detect leaks

Reagent Monitor Components

- RM01 monitor (100 to 2,000ml range)
- RM10 monitor (1,000 to 10,000ml range)
- RMH8 hub (allows extended cable lengths)
- Reagent Monitor API Kit
- Each Reagent Monitor includes Windows software & USB cable

Features

- Reports reagent level
- Low volume & over fill alarm
- Programmable alarm set points
- Scalable to facility-wide usage
- Easy-to-use software and API
- Outputs volume or weight via .CSV file

Compatible Labware

- Bottles / Carboys / Waste reservoirs